

A compact table top unit with high loading capacity and the performance of a big laboratory autoclave. The concept of LABOKLAV 25 is offering new technical solutions and a brand new design. The unit has a chamber lid made from temperature resistant safety glass, a transparent front cover made from safety glass, LED chamber illumination and coloured LED front lighting to visualize the cycle progress.

LABOKLAV 25 is extremely compact and space saving, but still offering the performance of a big laboratory autoclave:

- Sterilisation of liquids in open and slightly closed bottles
- Instrument sterilisation, also available with vacuum option offering pre-vacuum and drying vacuum
- Waste sterilisation at 134 °C

