

A compact table top unit with high loading capacity and the performance of a big laboratory autoclave. The concept of LABOKLAV 25 is offering new technical solutions and a brand new design. The unit has a chamber lid made from temperature resistant safety glass, a transparent front cover made from safety glass, LED chamber illumination and coloured LED front lighting to visualize the cycle progress.

LABOKLAV 25 is extremely compact and space saving, but still offering the performance of a big laboratory autoclave:

- Sterilisation of liquids in open and slightly closed bottles
- Instrument sterilisation, also available with vacuum option offering pre-vacuum and drying vacuum
- Waste sterilisation at 134 °C

LABOBLUE – air cooling, Version M

- air cooling
- for open vessels and slightly closed vessels

The unit is equipped with a high efficient re-cooling system to dramatically reduce cooling down times when sterilizing liquids. As a result the cooling down time can be reduced up to 50% (depending on loading volume and vessel size).

