

LABOARCHIEVE – data storage and recording made simple

- RS 485 integrated
- 4 MB Storage chip for storing process data integrated

The autoclaves can be equipped with different cycle recording systems. Data can be printed out with the integrated build-in printer (optional) or can be transferred onto a PC via the integrated interface (optional software required). Already the basic version of LABOKLAV has an integrated storage chip with 4 MB storage capacity. The cycle data is automatically stored and can be printed out or transferred to a PC at a later stage. A Minimum of 100 cycles are stored, more is possible but depends on the data volume per cycle.

