

LABOSELECT – one controller for all types

- foil key pad with clear and easy to understand symbols
- large LCD Display
- process data are displayed during programme run
- up to 20 pre-configured programmes
- code protected access level to programmes
- programmes can be individualized

A clear, easy to read and understand foil key pad eliminates operation errors and increases safety. Clear symbols guide through the menu, programmes can be opened fast. The structured big LCD display informs permanently during a programme cycle of all relevant process and programme parameters. In the unlikely event of an error a message is displayed.

Depending on the integrated options the 20 programmes are optimized for the machine configuration.

