

LABOVARIO – the universal installation concept

- space saving concept of installation, no matter whether installed as horizontal or vertical autoclave
- lid opening side is variable
- can be installed as horizontal or vertical autoclave, the unit remains the same

The series Laboklav has been designed to be very compact and space saving. All components are placed inside the housing. Additionally the unit offers various possibilities of installation. The lid can be opened to the right or to the left, the unit can be operated and loaded from the small side, it can be installed as a horizontal or vertical autoclave. All this without additional costs or any alterations. The compact design results in a lower loading height, making it easier to fill the autoclave with baskets and bottles when installed as vertical autoclave.

